



# Recombinant Mouse Angiopoietin-2 (Angpt2)

<b>Product Code</b>	CSB-YP001707MO
<b>Storage</b>	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
<b>Uniprot No.</b>	O35608
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Mus musculus (Mouse)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	<p>YS NFRKSV DSTG RRQYQVQNGP CSYTFLLPET DSCRSSSSPY            MSNAVQRDAP LDYDDSVQRL QVLENILENN TQWLMKLENY IQDNMCKEMV            EIQQNVVQNG TAVMIEIGTS LLNQTAQTR KLTDVEAQVL NQTTRLELQL            LQHSISTNKL EKQILDQTSE INKLQNKNSF LEQKVLDMEG KHSEQLQSMK            EQKDELQVLV SKQSSVIDEL EKKLV TATVN NSLLQKQQHD LMETVNSLLT            MMSSPNSKSS VAIRKEEQTT FRDCAEIFKS GLTTSGIYTL TFPNSTEEIK            AYCDMDVGGG GWTVIQHRED GSVDFQRTWK EYKEGFGSPL            GEYWLGNEFV SQLTGQHRVY LKIQLKDWEG NEAHSLYDHF YLAGEESNYR            IHLTGLTGTA GKISSISQPG SDFSTKSDSN DKCICKCSQM LSGGWWFDAC            GPSNLNGQYY PQQQNTNKNFN GIKWYWWKGS GYSLKATTMM IRPADF</p>
<b>Source</b>	Yeast
<b>Target Names</b>	Angpt2
<b>Protein Names</b>	Recommended name: Angiopoietin-2 Short name= ANG-2
<b>Expression Region</b>	19-496
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full Length of Mature Protein
<b>Target Details</b>	This protein is an antagonist of angiopoietin 1 (ANGPT1) and endothelial TEK tyrosine kinase (TIE-2, TEK). The encoded protein disrupts the vascular remodeling ability of ANGPT1 and may induce endothelial cell apoptosis. Three transcript variants encoding three different isoforms have been found for this gene.
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.