



Recombinant Arabidopsis thaliana 14-3-3-like protein GF14 epsilon (GRF10)

Product Code	CSB-BP343348DOA
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P48347
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	>85% (SDS-PAGE)
Sequence	MENEREKQVY LAKLSEQTER YDEMVEAMKK VAQLDVELTV EERNLVSVDGY KNVIGARRAS WRILSSIEQK EESKGNDENV KRLKNYRKRK EDELAKVCND ILSVIDKHLI PSSNAVESTV FFYKMKGDYY RYLAEFSSGA ERKEAADQSL EAYKAAVAAA ENGLAPTHPV RLGLALNFSV FYYEILNSPE SACQLAKQAF DDAIAELDSL NEESYKDSTL IMQLLRDNLTLWTSDLNEEG DERTKGADEP QDEN
Source	Baculovirus
Target Names	GRF10
Protein Names	Recommended name: 14-3-3-like protein GF14 epsilon Alternative name(s): General regulatory factor 10
Expression Region	1-254
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.