



Recombinant Probable S-adenosylmethionine synthase 4 (sams-4)

Product Code	CSB-YP343470CXY
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P50306
Product Type	Recombinant Protein
Immunogen Species	Caenorhabditis elegans
Purity	>85% (SDS-PAGE)
Sequence	MPQHKFLFTS ESVSEGHPDK MCDQISDAVL DAHLAQDPHA KVACETVTKT GMIMLCGEIT SKAVVDYQVL VRNVIKIGY DDSSKGFYK TCNVLVALEQ QSPEIAAGVH VDKSDDVGA GDQGIMFGYA TDETEEAMPL TLLLSHKLNR KLHELRRSGE LEWVRPDSKT QVTIEYASEG GACVPLRVHT VVISTQHSPD ISLDDLKEL IEKVIKAVIP ANLIDDKTIY HLNPCGSFII GGPMGDAGLT GRKIIVDTYG GWGAHGGGAF SGKDPTKVDR SAAYAARWVA KSLVKAGLCR RCLVQVSYAI GVAKPLSVMV FSFGTSALSE EDLLTIVNDN FDLRPGKIIK NLDLKRPIYE PTAENGHFGH NNFPEWQPRN LKISADMLAK SQGPAQPDVI GIAH
Source	Yeast
Target Names	sams-4
Protein Names	Recommended name: Probable S-adenosylmethionine synthase 4 Short name= AdoMet synthase 4 EC= 2.5.1.6 Alternative name(s): Methionine adenosyltransferase 4 Short name= MAT 4
Expression Region	1-404
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.