



Recombinant Arabidopsis thaliana Sedoheptulose-1,7-bisphosphatase, chloroplastic (At3g55800)

Product Code	CSB-EP343566DOA-B
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P46283
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	>85% (SDS-PAGE)
Sequence	A TKAKSNGAST VTKCEIGQSL EEFLAQATPD KGLRLLMCM GEALRTIAFK VRTASCGGTA CVNSFGDEQL AVDMLADKLL FEALQYSHVC KYACSEEVPE LQDMGGPVEG GFSVAFDPLD GSSIVDTNFT VGTIFGVWPG DKLTGITGGD QVAAAMGIYG PRTTYVLAVK GFPGTHEFLL LDEGKWQHVK ETTEIAEGKM FSPGNLRATF DNSEYSKLID YYVKEKYTLR YTGGMVPDVN QIIVKEKGIF TNVTSPTAKA KLRLLEFVAP LGLLIENAGG FSSDGHKSVL DKTIINLDDR TQVAYGSKNE IIRFEETLYG TSRLKNVPIG VTA
Source	E.coli
Target Names	At3g55800
Protein Names	Recommended name: Sedoheptulose-1,7-bisphosphatase, chloroplastic EC=3.1.3.37 Alternative name(s): SED(1,7)P2ase Sedoheptulose bisphosphatase Short name= SBPase
Expression Region	60-393
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.