



Recombinant *Oryza sativa* subsp. japonica 26S protease regulatory subunit 6A homolog (TBP1)

| | |
|--------------------------|--|
| Product Code | CSB-EP343582OFG-B |
| Storage | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |
| Uniprot No. | P46465 |
| Product Type | Recombinant Protein |
| Immunogen Species | <i>Oryza sativa</i> subsp. japonica (Rice) |
| Purity | >85% (SDS-PAGE) |
| Sequence | MSSPPAAAA AMAVDDADDD QLASMSTEDI VRATRLDNE TRVLKDELQR TNLEVESYKE KIKENQEKIK LNKQLPYLVG NIVEILEMNP EDEAEEDGAN IDLDSQRK GK CVVLKTSTRQ TIFLPVIGLV DPEKLPKPGDL VGVNKDSYLI LDTLPSEYDS RVKAMEVDEK PTEDYNDIGG LEKQIQELVE AIVLPMTHKD RFQKLGIRPP KGVLLYGPPG TGKTLMARAC AAQTNATFLK LAGPQLVQMF IGDGAKLVRD AFQLAKEKSP CIIFIDEIDA IGTKRFDSEV SGDREVQRTM LELLNQLDGF SSSDERIKVIA ATNRADILDP ALMRSGRLDR KIEFPHPSEE ARARILQIHS RKMNVNPDVN FEELARSTDD FNGAQLKAVC VEAGMLALRR DATEVTHEDF NEGIIQVQAK KKSSLNYYA |
| Source | E.coli |
| Target Names | TBP1 |
| Protein Names | Recommended name: 26S protease regulatory subunit 6A homolog Alternative name(s): Tat-binding protein homolog 1 Short name= TBP-1 |
| Expression Region | 1-429 |
| Notes | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week. |
| Tag Info | Tag type will be determined during the manufacturing process. |
| Protein Length | full length protein |
| Reconstitution | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference. |
| Shelf Life | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |