



Recombinant *Oryza sativa* subsp. japonica 26S protease regulatory subunit 6A homolog (TBP1)

Product Code	CSB-MP343582OFG
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P46465
Product Type	Recombinant Protein
Immunogen Species	<i>Oryza sativa</i> subsp. japonica (Rice)
Purity	>85% (SDS-PAGE)
Sequence	MSSPPAAAA AMAVDDADDD QLASMSTEDI VRATRLDNE TRVLKDELQR TNLEVESYKE KIKENQEKIK LNKQLPYLVG NIVEILEMNP EDEAEEDGAN IDLDSQRK GK CVVLKTSTRQ TIFLPVIGLV DPEKLPKPGDL VGVNKDSYLI LDTLPSEYDS RVKAMEVDEK PTEDYNDIGG LEKQIQELVE AIVLPMTHKD RFQKLGIRPP KGVLLYGPPG TGKTLMARAC AAQTNATFLK LAGPQLVQMF IGDGAKLVRD AFQLAKEKSP CIIFIDEIDA IGTKRFDSEV SGDREVQRTM LELLNQLDGF SSSDERIKVIA ATNRADILDP ALMRSGRLDR KIEFPHPSEE ARARILQIHS RKMNVNPDVN FEELARSTDD FNGAQLKAVC VEAGMLALRR DATEVTHEDF NEGIIQVQAK KKSSLNYYA
Source	Mammalian cell
Target Names	TBP1
Protein Names	Recommended name: 26S protease regulatory subunit 6A homolog Alternative name(s): Tat-binding protein homolog 1 Short name= TBP-1
Expression Region	1-429
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.