



Recombinant *Oryza sativa* subsp. *japonica* 60S ribosomal protein L16, mitochondrial (RPL16)

Product Code	CSB-BP343606OFG
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P46801
Product Type	Recombinant Protein
Immunogen Species	<i>Oryza sativa</i> subsp. <i>japonica</i> (Rice)
Purity	>85% (SDS-PAGE)
Sequence	MEKHLVMYLT RKSIMLLRKY PLVTEFQVSK CGSHIVKIRR DVLYPKRTKY SKYSKCRCSR GCEPDGTQLG FGRYGKTKSCR AGRLSYRAIE AARRATIGQF HRAMSGQFRR NCKIWVRVLA DLPITGKPAE VRMGRGKGNP TGWIARVSTG QIPFEMDGVS LSNARQAARL AAHKPCSSTK FVQWS
Source	Baculovirus
Target Names	RPL16
Protein Names	Recommended name: 60S ribosomal protein L16, mitochondrial
Expression Region	1-185
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.