



Recombinant Dermatophagoides farinae Mite allergen Der f 6 (DERF6)

Product Code	CSB-BP343783DCO
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P49276
Product Type	Recombinant Protein
Immunogen Species	Dermatophagoides farinae (American house dust mite)
Purity	>85% (SDS-PAGE)
Sequence	V VGGQDADLAE APFQISLLKD YLIMKRHMCG GSLISESTVV TAAHCTYGQK ASSLSVRYGT NQRTSSSYGD LKVKPIIQHE SYEQDQTQTD KTIILPNPV VPSTNVQMNE IETEDIVDGD KVTIYGWGLT DGNGKDLDPK LQKGSMTIVG NDRCNEKWGS INAIHPGMIC ALDKTQSGCN GDSGGPLVSA NRKLTGIVSW GPSKCPPGEY MSVFTRPKYY LDWITKNIV
Source	Baculovirus
Target Names	DERF6
Protein Names	Recommended name: Mite allergen Der f 6 EC= 3.4.21.- Alternative name(s): Allergen Der f VI DF5 Allergen= Der f 6
Expression Region	50-279
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.