



Recombinant Human RNA polymerase II elongation factor ELL2 (ELL2)

Product Code	CSB-MP007609HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	O00472
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MAAGGTGGLR EEQRYGLSCG RLGQDNITVL HVKLTETAIR ALETYQSHKN LIPFRPSIQF QGLHGLVKIP KNDPLNEVHN FNFYLSNVGK DNPQGSFDCI QQTFFSSGAS QLNCLGFIQD KITVCATNDS YQMTRETRMTQ AEEESRNRST KVIKPGGPYV GKRVRIRKAP QAVSDTVPER KRSTPMNPAN TIRKTHSSST ISQRPYRDRV IHLLALKAYK KPELLARLQK DGVNQKDKNS LGAILQQVAN LNSKDLSYTL KDYVFKELQR DWPGYSEIDR RSLESVLSRK LNPSQNAAGT SRSESPVCSS RDAVSSPQKR LLDSEFIDPL MNKKARISHL TNRPPTLNG HLNPTSEKSA AGLPLPPAAA AIPTPPPLPS TYLPISHPPQ IVNSNSNSPS TPEGRGTQDL PVDSFSQNDI IYEDQQDKYT SRTSLETLPV GSVLLKCPKP MEENHSMCHK KSKKKSCHKK EKDQIKKHI ETIEEKEEDL KREEEIAKLN NSSPNSSGGV KEDCTASMEP SAIELPDYLI KYIAIVSYEQ RQNYKDDFNA EYDEYRALHA RMETVARRFI KLDAQKRLS PGSKEYQNVH EEVLQEYQKI KQSSPNYHEE KYRCEYLHNK LAHIKRLIGE FDQQAESWS
Source	Mammalian cell
Target Names	ELL2
Protein Names	Recommended name: RNA polymerase II elongation factor ELL2
Expression Region	1-640
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.