



# Recombinant Solanum tuberosum Tubulin beta-2 chain (TUBB2)

<b>Product Code</b>	CSB-BP343963FIG
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	P46264
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Solanum tuberosum (Potato)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MREILHIQGG QCGNQIGSKF WEVICDEHGV DPTGRYKGT A AESDLQLERI NVYFNEASGG RYVPRAVLMD LEPGTMDSIR SGPYGQIFRP DNLVFGQSGA GNNWAKGHYT EGAELIDAVL DVVRKEAENC DCLQGFQVCH SLGGGTGSGM GTLLISKVRE EYPDRMMLTF SVFSPKVSD TVVEPYNATL SVHQLVENAD ECMVLDNEAL YDICFRTLKL TTPSFGDLNH LISATMSGVT CCLRFPGQLN SDLRKLAVNL IPLPRLHFFM VGFAPLTSRG SQQYISLTPV ELTQQMWD AK NMMCAADPRH GRYLTASAMF RGKMSTKEVD EQMINVQKN SSYFVEWIPN NVKSSVCDIP PTGLKMASTF VGNSTSIQEM FRRVSEQFTA MFRRKAFLHW YTGEGMDEME FTEAESNMND PVAEQYQQYQD ATADDEEEYD DEAADDHHQY ES
<b>Source</b>	Baculovirus
<b>Target Names</b>	TUBB2
<b>Protein Names</b>	Recommended name: Tubulin beta-2 chain Alternative name(s): Beta-2-tubulin
<b>Expression Region</b>	1-452
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.