



Recombinant *Saccharomyces cerevisiae* Cell wall mannoprotein CIS3 (CIS3)

Product Code	CSB-MP344023SVG
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P47001
Product Type	Recombinant Protein
Immunogen Species	<i>Saccharomyces cerevisiae</i> (strain ATCC 204508 / S288c) (Baker's yeast)
Purity	≥85% (SDS-PAGE)
Sequence	DVISQI GDGQVQATSA ATAQATDSQA QATTTATPTS SEKISSSASK TSTNATSSSC ATPSLKDSSC KNSGTLELTL KDGVLTDAGK RIGSIVANRQ FQFDGPPPQA GAIYAAGWSI TEDGYLALGD SDVIFYQCLSG NFYNLYDQNV AEQCSAIHLE AVSLVDC
Source	Mammalian cell
Target Names	CIS3
Protein Names	Recommended name: Cell wall mannoprotein CIS3 Alternative name(s): Covalently-linked cell wall protein 5/11 Protein with internal repeats 4 Soluble cell wall protein 8
Expression Region	65-227
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.