



Recombinant *Saccharomyces cerevisiae* Protein PRY1 (PRY1)

Product Code	CSB-YP344027SVG
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P47032
Product Type	Recombinant Protein
Immunogen Species	<i>Saccharomyces cerevisiae</i> (strain ATCC 204508 / S288c) (Baker's yeast)
Purity	>85% (SDS-PAGE)
Sequence	A PAVVTVTEHA HEAAVVTVQG VVYVENGQTR TTYETLAPAS TATPTSTATA LVAPPVAPSS ASSNSDVVLS ALKNLASVWG KTTDSTTTLT SSESTSQSLA QATTTSTPAA ASTTSTPAAT TTTSQAAATS SASSSDSDL SDFASSVLAEH NKKRALHKDT PALSWSDTLA SYAQDYADNY DCSGTLTHSG GPYGENLALG YDGPAAVDAW YNEISNYDFS NPGFSSNTGH FTQVVWKSTT QVGCIGKTCG GAWGDYVICS YDPAGNYEGE YADNVEPLA
Source	Yeast
Target Names	PRY1
Protein Names	Recommended name: Protein PRY1 Alternative name(s): Pathogen related in Sc 1
Expression Region	20-299
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.