



Recombinant Zea mays Phytoene dehydrogenase, chloroplastic/chromoplastic (PDS1)

Product Code	CSB-EP344176ZAX-B
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P49086
Product Type	Recombinant Protein
Immunogen Species	Zea mays (Maize)
Purity	>85% (SDS-PAGE)
Sequence	LQVV VAGAGLAGLS TAKYLADAGH KPILLEARDV LGGKVAAWKD EDGDWYETGL HIFFGAYPNI QNLFGELRIE DRLQWKEHSM IFAMPNKPGE FSRFDFPETL PAPINGIWAI LRNNEMLTWP EKVKFAIGLL PAMVGGQPYV EAQDGLTVSE WMKKQGV PDR VNDEVFIAMS KALNFINPDE LSMQCILIAL NRFLQEKHGS KMAFLDGNPP ERLCMPIVDH IRSRGGEVRL NSRIKKIELN PDGTVKHFAL SDGTQITGDA YVCATPVDIF KLLVPQEWSE ITYFKKLEKL VGVPVINVHI WFDRKLNNTY DHLLFSRSSL LSVYADMSVT CKEYYDPNRS MLELVFAPAD EWIGRSDTEI IDATMEELAK LFPDEIAADQ SKAKILKYHI VKTPRSVYKT VPNCEPCRPL QRSPIEGFYI AGDYTKQKYL ASMEGAVLSG KLCAQSIVQD YSRLALRSQK SLQSGEVPVP S
Source	E.coli
Target Names	PDS1
Protein Names	Recommended name: Phytoene dehydrogenase, chloroplastic/chromoplastic EC= 1.14.99.- Alternative name(s): Phytoene desaturase
Expression Region	97-571
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.