



# Recombinant *Saccharomyces cerevisiae* U4/U6 snRNA-associated-splicing factor PRP24 (PRP24)

<b>Product Code</b>	CSB-EP344233SVG-B
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	P49960
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	<i>Saccharomyces cerevisiae</i> (strain ATCC 204508 / S288c) (Baker's yeast)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MEYGHHARPD SKRPLDEGSP AAAGLTSKKA NEALTRNREL TTVLVKNLPK SYNQNKVYKY FKHC GPIIHV DVADSLKKNF RFARIEFARY DGALAAITKT HKVVGQNEII VSHLTECTLW MTNFPPSYTQ RNIRDLLQDI NVVALSIRLP SLRFNTSRRF AYIDVTSKED ARYCVEKLNK LKIEGYTLVT KVSNPLEKSK RTDSATLEGR EIMIRNLSTE LLDENLLRES FEGFGSIEKI NIPAGQKEHS FNNCCAFMVF ENKDSAERAL QMNRSLG NR EISVSLADKK PFLERNEVKR LLASRNSKEL ETLICLFPLS DKVSPSLICQ FLQEEIHINE KDIRKILLVS DFNGAIIIFR DSKFAAKMLM ILNGSQFQ GK VIRSGTINDM KRYYNQ QNH SMKHVKPSCI NMMEKGP NLQ VKKKIPDKQE QMSNDDFRKM FLGE
<b>Source</b>	E.coli
<b>Target Names</b>	PRP24
<b>Protein Names</b>	Recommended name: U4/U6 snRNA-associated-splicing factor PRP24 Short name= U4/U6 snRNP protein
<b>Expression Region</b>	1-444
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.