



Recombinant *Deinococcus radiodurans* L-lactate dehydrogenase (ldh)

Product Code	CSB-EP344303DJA-B
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P50933
Product Type	Recombinant Protein
Immunogen Species	<i>Deinococcus radiodurans</i> (strain ATCC 13939 / DSM 20539 / JCM 16871 / LMG 4051 / NBRC 15346 / NCIMB 9279 / R1 / VKM B-1422)
Purity	>85% (SDS-PAGE)
Sequence	MKVGVVGTGF VGSTAAFALV LRGSCSELVL VDRDEDRAQA EAEDIAHAAP VSHGTRVWHG GHSELADAQV VILTAGANQK PGESRLDLLE KNADIFRELV PQITRAAPDA VLLVTSNPVD LLTDLATQLA PGQPVIGSGT VLDSARFRHL MAQHAGVDGT HAHGYVLGEH GDSEVLAWSS AMVAGMPVAD FMQAQNLAWN EQVRAKIDEG TRNAAASIIIE GKRAYYYGIG AALARITEAV LRDRRAVLTV SAPTPEYGVV LSLPRVVGRQ GVLSTLHPKL TGDEQQKLEQ SAGVLRGFKQ QLGL
Source	E.coli
Target Names	ldh
Protein Names	Recommended name: L-lactate dehydrogenase Short name= L-LDH EC= 1.1.1.27
Expression Region	1-304
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.