



Recombinant Arabidopsis thaliana Homeobox-leucine zipper protein HAT2 (HAT2)

Product Code	CSB-MP344393DOA
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P46601
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	>85% (SDS-PAGE)
Sequence	MMMGKEDLGL SLSLGFSSQNH NPLQMNLNPN SSSLNNLQRL PWNQTFDPTS DLRKIDVNSF PSTVNCEEDT GVSSPNSTIS STISGKRSER EGISGTGVGS GDDHDEITPD RGYSRGTSDE EEDGGGETSRK KLRLSKDQSA FLEETFKEHN TLNPKQKLAL AKKLNLTARQ VEVWFQNRRA RTKLKQTEVD CEYLKRCVEK LTEENRRLQK EAMELRTLKL SPQFYGQMTP PTTLIMCPSC ERVGGPSSSN HHHNHRPVS I NPWVACAGQV AHGLNFEALR PRS
Source	Mammalian cell
Target Names	HAT2
Protein Names	Recommended name: Homeobox-leucine zipper protein HAT2 Alternative name(s): Homeodomain-leucine zipper protein HAT2 Short name= HD-ZIP protein 2
Expression Region	1-283
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.