



Recombinant Mouse Beta-tectorin (Tectb)

Product Code	CSB-YP023371MO
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	O08524
Product Type	Recombinant Protein
Immunogen Species	Mus musculus (Mouse)
Purity	≥85% (SDS-PAGE)
Sequence	KSC TPNKADVILV FCYPKTIITK IPECPYGWEV HQLALGGLCY NGVHEGGYYQ FVIPDLSPKN KSYCGTQSEY KPPIYHFYSH IVSNDSTVIV KNQPVNYSFS CTYHSTYLVN QAAFDQRVAT VHVKNGSMGT FESQLSLNFY TNAKFSTKKE APFVLETSEI GSDLFAGVEA KGLSVRFKVV LNCSWATPSA DFMYPLQWQL INKGCPDET VLVHENGKDH RATFQFNAFR FQNIPKLSKV WLHCETFICD SEKLSCPVNC DKRKRMLRDQ TGGVLVVELS LRSRA
Source	Yeast
Target Names	Tectb
Protein Names	Recommended name: Beta-tectorin
Expression Region	18-305
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.