



Recombinant *Saccharomyces cerevisiae*

Uncharacterized protein YJR154W (YJR154W)

Product Code	CSB-EP344442SVG
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P47181
Product Type	Recombinant Protein
Immunogen Species	<i>Saccharomyces cerevisiae</i> (strain ATCC 204508 / S288c) (Baker's yeast)
Purity	≥85% (SDS-PAGE)
Sequence	MNTDSHNLSE PYNIGGQKYI NMKKKEDLGV CQPGLTQKAF TVEDKFDYKA IIEKMEVYGL CVVKNFIETS RCDEILKEIE PHFYRYESWQ GSPFPKETTV ATRSVLHSST VLKDVVCDRM FCDISKHFLN EENYFAAGKV INKCTSDIQL NSGIVYKVGA GASDQGYHRE DIVHHTTHQA CERFQYGTET MVGLGVAFTD MNKENGSTRM IVGSHLWGPB DSCGNFDKRM EFHVNVAKGD AVLFLGSLYH AASANRTSQD RVAGYFFMTK SYLKPEENLH LGTDLRVFKG LPLEALQLLG LGISEPFCGH IDYKSPGHLI SSSLFENDIE KGY YGETIRV NYGSTQ
Source	E.coli
Target Names	YJR154W
Protein Names	Recommended name: Uncharacterized protein YJR154W
Expression Region	1-346
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.