



Recombinant Schizosaccharomyces pombe Zinc finger protein zfs1 (zfs1)

Product Code	CSB-MP344496SXV
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P47979
Product Type	Recombinant Protein
Immunogen Species	Schizosaccharomyces pombe (strain 972 / ATCC 24843) (Fission yeast)
Purity	>85% (SDS-PAGE)
Sequence	MVYSPMSRPQ VPLAFRQWPP YNYKSDPLVN SKLSQSTTSV NAGPSLISPS FLDSYANSSS LLHKPSTNLG SLKTSSLLAS DEVFPSSVMP QLRPLDSSLS VSPEMDGWPW HHNSSSNPQG YAWTPSLLSS NATSYLHSGS SPHGNTSNHP SPISSLES LP SRSSTGSGSL DFSGLANLHD DSKSLAMSLN MAGVPVSUDE NSQTSFPFVH GQPESTMSRK PKLCVQSKSM TNIRNSVAKP SLVRQSHSAG VIPIKPTASN ASIRNAPSNL SKQFSPSGNS PLTEASKPFV PQPSAAGDFR QAKGSASHPH GSGSSNGVAP NGKRALYKTE PCKNWQISGT CRYGSKCQFA HGNQELKEPP RHPKYKSERC RSFMMYGYCP YGLRCCFLHD ESNAQKSATI KQSP
Source	Mammalian cell
Target Names	zfs1
Protein Names	Recommended name: Zinc finger protein zfs1 Alternative name(s): Multicopy suppressor of overexpressed cyr1 protein 4
Expression Region	1-404
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.