



Recombinant Zea mays 14-3-3-like protein GF14-6 (GRF1)

Product Code	CSB-BP344569ZAX
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P49106
Product Type	Recombinant Protein
Immunogen Species	Zea mays (Maize)
Purity	>85% (SDS-PAGE)
Sequence	MASAELSREE NVYMAKLAEQ AERYEEMVEF MEKVAKTVDS EELTVEERNL LSVAYKNVIG ARRASWRIIS SIEQKEEGRG NEDRVTLIKD YRGKIETELT KICDGILKLL ETHLVPSSTA PESKVFYLMK KGDYYRYLAE FKTGAERKDA AENTMVAYKA AQDIALAELA PTHPIRLGLA LNFSVFYIEI LNTPDRACSL AKQAFDEAIS ELDTLSEESY KDSTLIMQLL RDNLTWLTSD ISEDPAAEIR EAPKRDSSEG Q
Source	Baculovirus
Target Names	GRF1
Protein Names	Recommended name: 14-3-3-like protein GF14-6
Expression Region	1-261
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.