



Recombinant Arabidopsis thaliana 60S ribosomal protein L4-2 (RPL4D)

Product Code	CSB-MP344608DOA
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P49691
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	>85% (SDS-PAGE)
Sequence	MVASAAARPL VTVQGLDGDM STDQSTTVTL PDVMTAPVRP DIVNFVHAQI SNNSRQPYAV SKKAGHQ TSA ESWG TGRAVS RIPRVPGGGT HRAGQA AFGN MCRGGRMFAP TKIWRRWHRR VNVNMKRHAI VSAIAATAVP ALVMARGHKI ENVPEMPLVV SDSAEAVEKT SAAIKVLKQI GAYDDAEKAK NSIGIRPGKG KMRNRRYISR KGPLVVFGE GAKIVKAFRN LPGVELCHVE RLNLLKLAPG GHLGRFVIWT KSAFEKLESI YGSFEKPSEK KKG YVLPRAK MVNADLARI NSDEVQSVVN PIKDGSKRAV LKKNPLKLN VMFKLN PYAK TAKRMSLLAE ASRVKAKKEK LEKKRKVVTK EEAQA IKAAG KAWYQTMISD SDYTEFDNFT KWLGASQ
Source	Mammalian cell
Target Names	RPL4D
Protein Names	Recommended name: 60S ribosomal protein L4-2 Alternative name(s): L1
Expression Region	1-407
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.