



# Recombinant *Vibrio harveyi* Chromosomal replication initiator protein DnaA (dnaA)

|                          |  |
|--------------------------|--|
| <b>Product Code</b>      | CSB-YP344622VFD  |
| <b>Storage</b>           | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.<br>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.   |
| <b>Uniprot No.</b>       | P49996   |
| <b>Product Type</b>      | Recombinant Protein  |
| <b>Immunogen Species</b> | <i>Vibrio harveyi</i> ( <i>Beneckea harveyi</i> )  |
| <b>Purity</b>            | >85% (SDS-PAGE)  |
| <b>Sequence</b>          | MSSSLWLQCL QQLQEELPAT EFSMWVRPLQ AELNDNTLTL FAPNRFVLDW<br>VRDKYLN SIN RLLQEYCGND VPSLRFVGS RRVAAPKPAP TRTPADVAEE<br>SSAPAQLQAR KPVHKTWDDD AQVIADINHR SNVNPCHKFN NFVEGKSNQL<br>GLAAARQVSD NPGAAYNPLF LYGGTGLGKT HLLHAVGNAI VDNNPNAKVV<br>YMHSERFVQD MVKALQNNAI EEFKRYRYSV DALLIDDIQF FANKERSQEE<br>FFHTFNALLE GNQQIILTSD RYPKEISGVE DRLKSRFGWG LTVAIEPPEL<br>ETRVAILMKK AEDHQIHLAD EVAFFIAKRL RSNVRELEGA LNRVIANANF<br>TGRPITIDFV REALRDLLAL QEKLVTIDNI QKTVAEYYKI KVADLLSKRR<br>SRSVARPRQL AMALAKELTN HSLPEIGDAF GGRDHHTTVLH ACRKIEQLRE<br>ESHDIKEDYS NLIRTLSS |
| <b>Source</b>            | Yeast  |
| <b>Target Names</b>      | dnaA   |
| <b>Protein Names</b>     | Recommended name: Chromosomal replication initiator protein DnaA   |
| <b>Expression Region</b> | 1-468  |
| <b>Notes</b>             | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.  |
| <b>Tag Info</b>          | Tag type will be determined during the manufacturing process.  |
| <b>Protein Length</b>    | full length protein  |
| <b>Reconstitution</b>    | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.  |
| <b>Shelf Life</b>        | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.<br>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.   |