



Recombinant Zea mays Sex determination protein tasselseed-2 (TS2)

Product Code	CSB-EP344632ZAX
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P50160
Product Type	Recombinant Protein
Immunogen Species	Zea mays (Maize)
Purity	>85% (SDS-PAGE)
Sequence	MHASLASYAA AAMPALDLRP EIAHAHQPVM SPSHHGWDGN GATAVPTPMP KRLDGKVAIV TGGARGIGEA IVRLFAGHGA RVVIADIDDA AGEALASALG PQVSFVRCDV SVEDDVRRRAV DWALSRHGGR LDVYCINNAGV LGRQTRAARS ILSFDAAEFD RVLRVNALGA ALGMKHAARA MAPRRAGSIV SVASVAAVLG GLGPHAYTAS KHAIVGLTKN AACELRAHGV RVNVCVSPFGV ATPMLINAWR QGHDDATADA DRDLDLDLDV TVPSDQEVEK MEEVVRGLAT LKGPTLRPRD IAEAVLFLAS DEARYISGHN LVVDGGVTTS RNLIGL
Source	E.coli
Target Names	TS2
Protein Names	Recommended name: Sex determination protein tasselseed-2
Expression Region	1-336
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.