



Recombinant Insertion element IS2A uncharacterized 48.2 kDa protein

Product Code	CSB-EP344688ENL
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P51026
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli
Purity	>85% (SDS-PAGE)
Sequence	MIVLILVFRL VIGEQMIDVL GPEKRRRRTT QEKIAIVQQS FEPGMTVSLV ARQHGVAASQ LFLWRKQYQE GSLTAVAAGE QVVPASELAA AMKQIKELQR LLGKKTME NE LLKEAVEYGR AKKWIAHAPL IARGWGVSLV SRCLRV SRAQ LHVILRRRTDD WMDGRRSRHT DTDVLLRIH HVIGELPTYG YRRVWALLRR QAELDGMPAI NAKRVYRIMR QNALLERKP AVPPSKRAHT GRVAVKESNQ RWCSDFEFC CDNGERLRVT FALDCCDREA LHWAVTTGGF NSETVQDVML GAVERFRGND LPSSPVEWLT DNGSCYRANE TRQFARMLGL EPKNTAVRSP ESNGLAESFV KTIKRDIYISI MPKPDGLTAA KNLAEAFEHY NEWHPHSALG YRSPREYLRH GACNGLSDNR CLEI
Source	E.coli
Protein Names	Recommended name: Insertion element IS2A uncharacterized 48.2 kDa protein
Expression Region	1-424
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.