



Recombinant Arabidopsis thaliana Serine/threonine-protein kinase AFC1 (AFC1)

Product Code	CSB-BP344726DOA
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P51566
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	>85% (SDS-PAGE)
Sequence	MQSSVYRDKA SSIAMILETQ RNVEFPHRIV DKRPRKRPRL TWDAAPPLLP PPPPPTVFQP PLYYGPEFAS GLVPNFVYPN MFYNGLPRQG SPPWRPDDKD GHYVFFVVGDT LTPRYQILSK MGEFTFGQVL ECFDNKNKEV VAIKVIRSIN KYREAAMIEI DVLQRLTRHD VGGSRCVQIR NWFYRDNHIC IVFEKLGPSL YDFLRKNSYR SFPIDLVREL GRQLLESVAY MHDRLRIHTD LKPENILLVS SEYIKIPDYK FLSRPTKDGS YFKNLPKSSA IKLIDFGSTT FEHQDHNIV STRHYRAPEV ILGVGWNYPC DLWSIGCILV ELCSGEALFQ THENLEHLAM MERVLGPLPP HMVLRADRRS EKYFRRGAKL DWPEGATSRD SLKAVWKLPR LPNLIMQHVD HSAGDLIDLL QGLLRYPTE RFKAREALNH PFFTRSREQS IPPFNPNPFP FLYNQKN
Source	Baculovirus
Target Names	AFC1
Protein Names	Recommended name: Serine/threonine-protein kinase AFC1 EC= 2.7.12.1
Expression Region	1-467
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.