



Recombinant Danio rerio Transcriptional regulator Myc-A (myca)

Product Code	CSB-YP344765DIL
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P52160
Product Type	Recombinant Protein
Immunogen Species	Danio rerio (Zebrafish) (Brachydanio rerio)
Purity	>85% (SDS-PAGE)
Sequence	MPVSASLACK NYDYDYDSIQ PYFYFDNDDE DFYHHQQGQT QPSAPSEDIW KKFELLPTPP LSPSRRQSL S TAEQLEMVSE FLGDDVVSQS FICDDADYSQ SFIKSIIIQD CMWSGFSA A KLEKVVSERL ASLHAERKEL MSDSNSNRLN ASYLQDLSTS ASECIDPSV V FPYPLTECGK AGKVASQPM LVLDTPPNSS SSSGSDSEDE EEEDEEEEE EEEEEEEEE E EIDVVTVEK RQKRHETDAS ESRYSPPLVL KRCHVSTHQH NYAAHPSTRH DQPAVKRLRL EASNHSINS SSSNRHVKQR KCASPRTSDS EDNDKRRRTHN VLERQRRNEL KLSFFALRDE IPEVANNEKA AKVVILKKAT ECIHSMQLDE QRLLSIKEQL RRKSEQLKHR LQQLRSSH
Source	Yeast
Target Names	myca
Protein Names	Recommended name: Transcriptional regulator Myc-A Short name= c-Myc-A Short name= zc-Myc
Expression Region	1-408
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.