



Recombinant Solanum tuberosum Isoflavone reductase homolog

Product Code	CSB-EP344802FIG-B
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P52578
Product Type	Recombinant Protein
Immunogen Species	Solanum tuberosum (Potato)
Purity	>85% (SDS-PAGE)
Sequence	MAGKSKILFI GGTGYIGKFI VEASAKAGHD TFVLVRESTL SNPTKTKLID TFKSGVTFV HGDLYDHESL VKAIKQVDVV ISTVGHALLA DQVKLIAAIK EAGNVKRFFP SEFGNDVDRV HAVEPAKAAF NTKAQIRRVV EAEGIPFTYV ATFFFAGYSL PNLAQPGAAG PPNDKVVILG HGNTKAVFNK EEDIGTYTIN AVDDPKTLNK ILYIKPPHNI ITLNELVSLW EKKTGKNLER LYVPEEQVLK NIQEASVPMN VGLSIYHTAF VKGDHTNFEI EPSFGVEASE VYPDVKYTP DEILNQYV
Source	E.coli
Protein Names	Recommended name: Isoflavone reductase homolog EC= 1.3.1.- Alternative name(s): CP100
Expression Region	1-308
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.