



Recombinant Uncharacterized serine carboxypeptidase F13S12.6 (F13D12.6)

Product Code	CSB-EP344819CXY-B
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P52715
Product Type	Recombinant Protein
Immunogen Species	Caenorhabditis elegans
Purity	>85% (SDS-PAGE)
Sequence	ALITNLPGA PISNFKQYSG YYNVGTTKKNH MLHYWFVESQ SNPSTDPVLL WLTGGPGCSG LSALLTEWGP WNVNTDGATL RTNPYSWNKN ASILTLEAPA GVGYSYATDN NIATGDDQTA SENWEALVAF FNEFPQYKGN DFYVTGESYG GIYVPTLVQT ILDRQSQSHI NIKGLAINGN CVSANEGVDS LVNFLYHHGV VDQAKWEHMK TSCCHNDTDA CPWHSFSEFS ACGEFVEATQ QTAWNGGLNP YNMYADCIST SASFRFGMEY ERRFNKKYTP EVLGTVPCLD ESPTNYLNR QDVRKALGIP SSLPAWSICS NAISYGYKRQ YGDMTSRVLN AVNNNNLKMM LYNGDVDLAC NALMGQRFTD KLGLTSLSKK THFTVKGQIG GYVTQYKGSQ VTFATVRGAG HVMPTDKPAV AEHIIQSFLF NKAFL
Source	E.coli
Target Names	F13D12.6
Protein Names	Recommended name: Uncharacterized serine carboxypeptidase F13S12.6 EC=3.4.16.-
Expression Region	22-454
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.