



Recombinant *Saccharomyces cerevisiae* Cytosolic Fe-S cluster assembly factor NBP35 (NBP35)

Product Code	CSB-MP344831SVG
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P52920
Product Type	Recombinant Protein
Immunogen Species	<i>Saccharomyces cerevisiae</i> (strain ATCC 204508 / S288c) (Baker's yeast)
Purity	>85% (SDS-PAGE)
Sequence	MTEILPHVND EVLPAEYELN QPEPEHCPGP ESDMAGKSDA CGGCANKEIC ESLPKGPDPD IPLITDNLSG IEHKILVLSG KGGVGKSTFA AMLSWALSAD EDLQVGAMD L DICGPSLPHM LGCIKETVHE SNSGWTPVYV TDNLATMSIQ YMLPEDDSAI IWRGSKKNLL IKKFLKDVDW DKLDYLVIDT PPGTSDEHIS INKYMRESGI DGALVVTPQ EVALLDVRKE IDFCKKAGIN ILGLVENMSG FVCPNCKGES QIFKATTGGG EALCKELGIK FLGSVPLDPR IGKSCDMGES FLDNYPDSPA SSAVLNVVEA LRDAVGDV
Source	Mammalian cell
Target Names	NBP35
Protein Names	Recommended name: Cytosolic Fe-S cluster assembly factor NBP35 Alternative name(s): Nucleotide-binding protein 35
Expression Region	1-328
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.