



# Recombinant *Saccharomyces cerevisiae* Putative cysteine synthase (YGR012W)

<b>Product Code</b>	CSB-EP344862SVG-B
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	P53206
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	<i>Saccharomyces cerevisiae</i> (strain ATCC 204508 / S288c) (Baker's yeast)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	MSCSQNKTSV SLAWRECISI ASVLIGAYAS YKYYKLFKTR DIPRPKEGVE ELIGNTPLVK IRSCLKATGV NIYAKLELCN PAGSAKDRVA LNIKTAEEL GELVRGEPGW VFEGTSGSTG ISIAVVCNAL GYRAHISLPD DTSLEKLALL ESLGATVNVKVPASIVDPNQ YVNAACKACN ELKKSNGNIR AVFADQFENE ANWKVHYQTT GPEIAHQTKG NIDAFIAGCG TGGTITGVAK FLKERAKIPC HVVLADPQGS GFYNRVNYGV MYDYVEKEGT RRRHQVDIV EGIGLNRITH NFHMGEKFID ESIRVNDNQA IRMAKYL SVN DGLFVGSSTA INAVAAIQVA KTLPHGSNIV IACDSGSRH LSKFWKEAKE IDHDSVLEEV INI
<b>Source</b>	E.coli
<b>Target Names</b>	YGR012W
<b>Protein Names</b>	Recommended name: Putative cysteine synthase Short name= CSase EC= 2.5.1.47 Alternative name(s): O-acetylserine (thiol)-lyase Short name= OAS-TL O-acetylserine sulfhydrylase
<b>Expression Region</b>	1-393
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.