



Recombinant *Saccharomyces cerevisiae* Probable glycosidase CRH1 (CRH1)

Product Code	CSB-BP344875SVG
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P53301
Product Type	Recombinant Protein
Immunogen Species	<i>Saccharomyces cerevisiae</i> (strain ATCC 204508 / S288c) (Baker's yeast)
Purity	>85% (SDS-PAGE)
Sequence	ESTATADS TTAASSTASC NPLKTTGCTP DTALATSFSE DFSSSSKWFT DLKHAGEIKY GSDGLSMTLA KRYDNPSLKS NFYIMYGKLE VILKAANGTG IVSSFYLQSD DLDEIDIEWV GGDNTQFQSN FFSKGDTTTTY DRGEFHGVDT PTDKFHNYTL DWAMDKTTWY LDGESVRVLS NTSSEGY PQS PMYLMMGIWA GGDPD NAAGT IEWAGGETNY NDAPFTMYIE KVIVTDYSTG KKYTYGDQSG SWESIEADGG SIYGRYDQAQ EDFAVLANGG SSSSSTSSS TVSSASSTV SSSVSSTVSS SASSTVSSSV SSTVSSSSSV SSSSSTSPSS STATSSKTLA SSSVTTSSSI SSFEKQSSSS SKKTVASSST SESIISSTKT PATVSSTTRS TVAPTTQSS VSSDSPVQDK GGVATSSNDV TSSTTQISSK YTSTIQSSSS EASSTNSVQI SN
Source	Baculovirus
Target Names	CRH1
Protein Names	Recommended name: Probable glycosidase CRH1 EC= 3.2.-.- Alternative name(s): Congo red hypersensitive protein 1
Expression Region	23-482
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.