



Recombinant *Saccharomyces cerevisiae* Maltose fermentation regulatory protein MAL13 (MAL13)

Product Code	CSB-YP344880SVG
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P53338
Product Type	Recombinant Protein
Immunogen Species	<i>Saccharomyces cerevisiae</i> (strain ATCC 204508 / S288c) (Baker's yeast)
Purity	>85% (SDS-PAGE)
Sequence	MTLTKQTCAK QACDCCRIRR VKCDGKRPCS SCLQNSLDCT YLQPSRKRGPKSIRLRSLKR IAEVQRESGP NTIATAPVIY KRVPKKLIDQ CLRLYHDNLYVIWPLLSYDD LHKLLEEKYN DNYVYWFLTA LSAATLSDLQ TEIKSEEEVFTGKQLSNLC ISSCQQFDDL DNSNIFNIMT YYCLHRFAQ ISNARTSYRLCCEAVGLITV AGLHREETYG SLTFEEQQLR RKLYYLLMT ERYyaiylhcATSLDATIAP PQLELVTDPO LSMDSFLEMI RVFTVPGKCF FDALAADSTDASCTEESLKK IWNELHTTSS EIEPWSNGYI DISFSRHWIR ILAWKLAYQMRGSNFSLNAN NGQIPIEIAR DMLIDTYLTP ENLYDVHGPG VPKTLEIATALVDIVGQYD HNMKLEAWNVLHDVCKFAFS LNHYNNDMLK RFSTKQCQNALITLPIISKPLQ LDGYPKDNED IDP
Source	Yeast
Target Names	MAL13
Protein Names	Recommended name: Maltose fermentation regulatory protein MAL13
Expression Region	1-473
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.