



Recombinant Human Ubiquitin domain-containing protein UBFD1 (UBFD1)

Product Code	CSB-MP025495HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	O14562
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MAAAGAPDGM EEPGMDTEAE TVATEAPARP VNCLEAEAAA GAAAEDSGAA RGS LQPAPAQ PPGDPAAQAS VSNGEDAGGG AGRELVDLKI IWNKTKHDVK FPLDSTGSEL KQKIHSITGL PPAMQKVMYK GLVPEDKTLR EIKVTSGAKI MVGSTINDV LAVNTPKDAA QQDAKAEENK KEPLCRQKQH RKVLDKGKPE DVMP SVKGAQ ERLPTVPLSG MYNKS GGKVR LTFKLEQDQL WIGTKERTEK LPMGSIKNVV SEPIEGHEDY HMMAFQLGPT EASY YWVYVW PTQYVDAIKD TVLGKWQYF
Source	Mammalian cell
Target Names	UBFD1
Protein Names	Recommended name: Ubiquitin domain-containing protein UBFD1 Alternative name(s): Ubiquitin-binding protein homolog
Expression Region	1-309
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.