



Recombinant *Saccharomyces cerevisiae* Lipase 2 (TGL2)

Product Code	CSB-EP345017SVG
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P54857
Product Type	Recombinant Protein
Immunogen Species	<i>Saccharomyces cerevisiae</i> (strain ATCC 204508 / S288c) (Baker's yeast)
Purity	≥85% (SDS-PAGE)
Sequence	MKNDNKANDI IIDSVKVPDS YKPPKNPIVF CHGLSGFDKL ILIPSVFHLT NLISNSIVHN MAENFMQDDE DKSDNKYTNL LEIEYWIGVK KFLQSKGCTV ITTKVPGFGS IEERAMALDA QLQKEVKKIE SKDKRHSLNL IAHSMSGGLDC RYLICNIKNR NYDILSLTTI STPHRGSEMA DYVVDLFENL NALRVSQKIL PICFYQLTTA YMKYFNLVTP NSPKVSYFSY GCSFVPKWYN VFCTPWKIVY ERSKGCPCNDG LVTINSSKWG EYRGTLDKMD HLDVINWKNK LQDDWSKFFR TTTVGEKVDI LNFYLRKITDD LARKGF
Source	E.coli
Target Names	TGL2
Protein Names	Recommended name: Lipase 2 EC= 3.1.1.3 Alternative name(s): Triacylglycerol lipase 2
Expression Region	1-326
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.