



# Recombinant Arabidopsis thaliana IAA-amino acid hydrolase ILR1 (ILR1)

<b>Product Code</b>	CSB-EP345031DOA-B
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	P54968
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Arabidopsis thaliana (Mouse-ear cress)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	GSYDSGS GLESLARGML HSAKDPEFFE WMRGIRRKIH ENPETGFQEF KTSQLVRDEL DSLGVKYKYP VAKTGVVAWI GSCSKPVFGL RADMDALPLQ ELVEWESKSK VDGKMHACGH DTHVAMLLGA AKLLQTTKHL IKGTVKLVFQ PGEEGYAGAY EMLKDEILDD LDGILSVHVF PSIPSGGIGS RPGTVLAGAG LFTVTVHGQG SHAATPHFSK DPVLAASSAV VALQQIVSRE LDPLEAGVVT VGYIEGGHAQ NVIPQSAKFG GTFRSLSDNG LLFIQRRIKE ISEAQASVYR CKAEVNFEEK KPSLHPVMNN DEGLYEHGKK VAEAMIGKNN FHDFPVTMGG EDFSFFTQKT KAAIFVLGVK NETLGAGKPL HSPYFFVDEE ALPVGAALHA AMAVSYLDEH GHSHEEEVKS EL
<b>Source</b>	E.coli
<b>Target Names</b>	ILR1
<b>Protein Names</b>	Recommended name: IAA-amino acid hydrolase ILR1 EC= 3.5.1.-
<b>Expression Region</b>	24-442
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full Length of Mature Protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.