



Recombinant Uncharacterized serine carboxypeptidase F41C3.5 (F41C3.5)

Product Code	CSB-MP345228CXY
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P52717
Product Type	Recombinant Protein
Immunogen Species	Caenorhabditis elegans
Purity	≥85% (SDS-PAGE)
Sequence	EEIK DLPGLDFEPN FKHYSGFFQV SDNHVLHYWF VESQNEPSND PLIFWFNGGP GCSSLDGLLN EMGPYVANED GKTLRENEYS WNKMASVVYI ESPAGVGYSY ATDGNITTND DLTSLENYEA VKQFFTEFPQ FRHHQTFIMG ESYGGVYVPT LTARIVDGQK DFPINLKGMA LGNGYVNEKL NIDTSVRFAY GHGLIDEKIW NTLERDCCRG CIDSCDLTQV TGHCATLVED IFQFLWFGGL NPYDLYRDCD PNPSVNSKRM KHMLRGVAPA MAHFDELLKN QTKTSLYQFL KNKSQSQKPL KADVPCLNLT EMLSYMNNPK VRKAIHIPFN LGKWDICSDK VTTTTYQKQYT DMTPFIKKIV KNHVRVLLYY GDTDMACNFM MGQQFSDQLG LRRTLKKTWP KYDRQIAGFK TLFDGLSFIT IRGAGHMAPQ WRAPQMYYAV QQFLNNHPL
Source	Mammalian cell
Target Names	F41C3.5
Protein Names	Recommended name: Uncharacterized serine carboxypeptidase F41C3.5 EC=3.4.16.-
Expression Region	17-469
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.