



Recombinant Arabidopsis thaliana Pyruvate dehydrogenase E1 component subunit alpha-1, mitochondrial (At1g59900)

Product Code	CSB-BP345237DOA
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P52901
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	≥85% (SDS-PAGE)
Sequence	TIETSLPF TAHLCDPPSR SVESSSQELL DFFRTMALMR RMEIAADSLY KAKLIRGFCH LYDGQEAVAI GMEAAITKGD AITAYRDHC IFLGRGGS LH EVFSELMGRQ AGCSKGKGGG MHFYKKESS YGGHGIVGAQ VPLGCGIAFA QKYNKEEAVT FALYGDGAAN QGQLFEALNI SALWDLPAIL VCENNHYGMG TAEWRAAKSP SYYKRGDYVP GLKVDGMDAF AVKQACKFAK QHALEKGP II LEMDTYRYHG HSMSDPGSTY RTRDEISGVR QERDPIERIK KLVLSHDLAT EKELKDMEKE IRKEVDDAIA KAKDCPMPEP SELFTNVYVK GFGTESFGPD RKEVKASLP
Source	Baculovirus
Target Names	E1
Protein Names	Recommended name: Pyruvate dehydrogenase E1 component subunit alpha-1, mitochondrial Short name= PDHE1-A EC= 1.2.4.1
Expression Region	33-389
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.