



Recombinant *Saccharomyces cerevisiae*

Uncharacterized protein YGR111W (YGR111W)

Product Code	CSB-BP345273SVG
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P53265
Product Type	Recombinant Protein
Immunogen Species	<i>Saccharomyces cerevisiae</i> (strain ATCC 204508 / S288c) (Baker's yeast)
Purity	>85% (SDS-PAGE)
Sequence	MTASSNDDDL IFECYSDEL KRWTHLANAK AWKGILTVQQ YADREQLLGS SEISQKNKSN EMMTKPKSY QWLGQKYFVL KDRSLPDNGK FSQVVSSCET LNRIGYCIHP GSNGKIEPAL IVCIGGVFTF ENHRGKGYAK KMIIKLNEFY DKIRDDANTV LELKNLVINL YSEVGEYYSALGYESMHVPL HRISKLDEL ERYCGEDDDH DGKYLGFDDY RGLVGLHETQ FKESLLSLHK ENPEKVFVTV APDFDIFTWF QYRDLFIMNK SGRKAQQNLF FGYALSDNSH IIVHHNWNGD SLIIVKIHIP EETFQRKELK LKLLRKAIE ETKLHGLQEL EFWDEEPIK KYPQLFQLLT ELENESKVFS ENGSISAVRP PKGYTAEQVI WDNNTKFCWF
Source	Baculovirus
Target Names	YGR111W
Protein Names	Recommended name: Uncharacterized protein YGR111W
Expression Region	1-400
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.