




# Recombinant Human Mucin-2 (MUC2), partial

<b>Product Code</b>	CSB-YP015222HU
<b>Relevance</b>	Coats the epithelia of the intestines, airways, and other mucus membrane-containing organs. Thought to provide a protective, lubricating barrier against particles and infectious agents at mucosal surfaces. Major constituent of both the inner and outer mucus layers of the colon and may play a role in excluding bacteria from the inner mucus layer.
<b>Abbreviation</b>	Recombinant Human MUC2 protein, partial
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q02817
<b>Alias</b>	Intestinal mucin-2
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	≥ 90% as determined by SDS-PAGE.
<b>Sequence</b>	VCSTWGNFHYKTFDGDVFRFPGLCDYNFASDCRGSYKEFAVHLKRGPGQAE APAGVESILLTIKDDTIYLTRHLAVLNGAVVSTPHYSPGLLIEKSDAYTKVYSRA GLTLMWNREDALMLELDTKFRNHTCGLCGDYNGLQSYSEFLSDGVLFSPLEF GNMQKINQPDVVCEDPEEEVAPASCSEHRAECERLLTAEAFADCQDL
<b>Research Area</b>	Signal Transduction
<b>Source</b>	Yeast
<b>Target Names</b>	MUC2
<b>Protein Names</b>	Recommended name: Mucin-2 Short name= MUC-2 Alternative name(s): Intestinal mucin-2
<b>Expression Region</b>	36-240aa
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	N-terminal 6xHis-tagged
<b>Mol. Weight</b>	24.8kDa
<b>Protein Length</b>	Partial
<b>Image</b>	



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

### Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

### Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.