



Recombinant *Drosophila melanogaster* Putative peptidyl-prolyl cis-trans isomerase dodo (dod)

Product Code	CSB-EP345377DLU
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P54353
Product Type	Recombinant Protein
Immunogen Species	<i>Drosophila melanogaster</i> (Fruit fly)
Purity	>85% (SDS-PAGE)
Sequence	MPDAEQLPDG WEKRTSRSTG MSYYLNMYTK ESQWDQPTPEP AKKAGGGSAG GGDAPDEVHC LHLLVKHKGS RRPSSWREAN ITRTKEEAQL LLEVYRNKIV QQEATFDELA RSYSDCSSAK RGGDLGKFGR GQMQAAFEDA AFKLNVNQLS GIVDSDSGLH IILRKA
Source	E.coli
Target Names	dod
Protein Names	Recommended name: Putative peptidyl-prolyl cis-trans isomerase dodo Short name= PPIase dodo EC= 5.2.1.8 Alternative name(s): Rotamase dodo
Expression Region	1-166
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.