



# Recombinant Human Eukaryotic translation initiation factor 3 subunit D (EIF3D)

<b>Product Code</b>	CSB-EP007533HU-B
<b>Storage</b>	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
<b>Uniprot No.</b>	O15371
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	<p>MAKFMTPIVIQ DNPSGWGPCA VPEQFRDMPY QPFSKGDRLG  KVADWTGATY QDKRYTNKYS SQFGGGSQYA YFHEEDESSF  QLVDTARTQK TAYQRNRMRF AQRNLRRDKD RRNMLQFNLQ ILPKSAKQKE  RERIRLQKKF QKQFGVRQKW DQKSQKPRDS SVEVRSDWEV  KEEMDFPQLM KMRYLEVSEP QDIECCGALE YYDKAFDRIT TRSEKPLRSI  KRIFHTVTTT DDPVIRKLAK TQGNVFATDA ILATLMSCTR SVYSWDIVVQ  RVGSKLFFDK RDNSDFDLLT VSETANEPQ DEGNSFNPR NLAMEATYIN  HNFSQQCLRM GKERYNFPNP NPFVEDDMDK NEIASVAYRY RRWKLGDID  LIVRCEHDGV MTGANGEVSF INIKTLNEWDRHRCNGVDWR QKLDSSQRGAV  IATELKNNSY KLARWTCCAL LAGSEYLKLG YVSRVHVKDS SRHVILGTQQ  FKPNEFASQI NLSVENAWGI LRCVIDICMK LEEGKYLIK DPNKQVIRVY  SLPDGTFSSD EEEEEEEEEEE EEEEEEEET</p>
<b>Source</b>	E.coli
<b>Target Names</b>	EIF3D
<b>Protein Names</b>	Recommended name: Eukaryotic translation initiation factor 3 subunit D Short name= eIF3d Alternative name(s): Eukaryotic translation initiation factor 3 subunit 7 eIF-3-zeta eIF3 p66
<b>Expression Region</b>	1-548
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full length protein
<b>Target Details</b>	Eukaryotic translation initiation factor-3 (eIF3), the largest of the eIFs, is a multiprotein complex composed of at least ten nonidentical subunits. The complex binds to the 40S ribosome and helps maintain the 40S and 60S ribosomal subunits in a dissociated state. It is also thought to play a role in the formation of the 40S initiation complex by interacting with the ternary complex of eIF2/GTP/methionyl-tRNA, and by promoting mRNA binding. This protein is the major RNA binding subunit of the eIF3 complex.
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a



concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

### Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.