



# Recombinant Peptidyl-prolyl cis-trans isomerase 8 (cyn-8)

<b>Product Code</b>	CSB-EP345563CXY-B
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	P52016
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Caenorhabditis elegans
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	MPPEVRGNKR AFFDISINGE PAGRIVFSLW NHCCPRTVEN FRAFCTGELG KMNGHYASYQ GSVFHRVIKG FMIQGGDITH GNGTGGYSIY GRTFDDENLA LKHKKPYLLS MANRGPDTNG SQFFITSEEV PHLDGKHCVF GEVIKGV EKV KAIENLETGN EDKPVCKVEI THCGEMVRKG DVVGNAGAAA SEPEPTGVKM EDEPKTRNWL MRYSKSPESR EEKKKDKHGR EEKRDRRRRS NDRHGRDRRS RRSRSQRDRN RRRDDRSGRD GRVGRNERDD RSGRDERRDD RRSSAPVNVK KKTQMKSPNF VHAPTKDGIK VRGRGKLTFL GGRTRSTTPP HWKREESKKL TLEAHQKRQE ELEERQKRAA QREKEEAERN ARLEKERVER QNQRELERQK QEAERELKRQ QENEESGNRR RRRSSRSRGSE NDESELEEPK NRSASPRKRS SSSSSSSSSG HSSDSD
<b>Source</b>	E.coli
<b>Target Names</b>	cyn-8
<b>Protein Names</b>	Recommended name: Peptidyl-prolyl cis-trans isomerase 8 Short name= PPIase 8 EC= 5.2.1.8 Alternative name(s): Cyclophilin-8 Rotamase 8
<b>Expression Region</b>	1-466
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.