



Recombinant Escherichia coli Uncharacterized protein yfjK (yfjK)

Product Code	CSB-BP345570ENV
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P52126
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli (strain K12)
Purity	≥85% (SDS-PAGE)
Sequence	MTEIYEQAKH SLQGEDFSSF NYLFAVNKLL SNPVSYDLGR DLIVRALDSR ERFSEHTTIL KNMVRKSGLF PYLKKEFTSL TPDDLRLVLEL YRTPFSDGYV FHSMQFHIFD LLKSGQNVVL SAPTSMGKSA IVDSLLGMGT LKRLVLVPT VALADETRRR LQERFGDRYQ IIHSSQVCH SDQAVYVLTQ ERVNERDDIV DIDLFFVIDEF YKLAFRQLKS GDIDHQDERV IELNIALSKL LKVSQRQFYLT GPFVNSIRGL EKLGYPTFTV STDFNTVALD VKTFGIKAND DKAKLKALGE IAHACVDATI IYCKSPTVAG LVARELIRLG HGTPTENPHV DWVSEEFDAD WDYTVLRNG IGLHFGALPR ALQQYTADQF NAGKLRFLLC TSTIIEGVNT IAKNVVIYDN RDGTRSIDKF THGNIKGRAG RMGVHVFVGKI FCLEEIPEDN LNQEVDIPLG IQGIDTPINL LASVQPDHLS EFSQDRFDEV FINDRVSIDL VKKHSYFRVE QFEMLQSMFE MMDDNEFSSL VFHWTPATNF LKTFAKIAR LVPHTFSRNG VPVKPTDVM IAKLAGYLSAE SYSEYLKNQI DYARQWISEG EKRTLSIALN NDLKLITNTF GYTLPKVLSL MEDVVKHHAV KRGIRSKVDY THVKLAFESF HLPPGVNALE EIGIPIQTLH RLVDLLEFSD EADVDELSQY LRDTQDIWSR SIGYVDQMF IRRALGIRRH
Source	Baculovirus
Target Names	yfjK
Protein Names	Recommended name: Uncharacterized protein yfjK
Expression Region	1-729
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.