



Recombinant Zea mays Isoflavone reductase homolog IRL (IRL)

Product Code	CSB-MP345612ZAX
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P52580
Product Type	Recombinant Protein
Immunogen Species	Zea mays (Maize)
Purity	≥85% (SDS-PAGE)
Sequence	MASEKSKILV VGGTGYLGRH VVAASARLGH PTSALVRDTA PSDPAKAALL KSFQDAGVTL LKGDLYDQAS LVSAVKGADV VISVLGSMQI ADQSRLVDAI KEAGNVKRFF PSEFGLDVDR TGIVEPAKSI LGAKVGIRRA TEAAGIPYTY AVAGFFAGFG LPKVGQVLAP GPPADKAVVL GDGDTKAVFV EEGDIATYTV LAADDPRAEN KVLVIKPPAN TLSHNELLSL WEKKTGKTFR REYVPEEAVL KQIQESPIPL NIILAIGHAA FVRGEQTGFE IDPAKGVDAS ELYPDVKYTT VDEYLNRFLL
Source	Mammalian cell
Target Names	IRL
Protein Names	Recommended name: Isoflavone reductase homolog IRL EC= 1.3.1.-
Expression Region	1-309
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.