



Recombinant *Saccharomyces cerevisiae* J-type co-chaperone JAC1, mitochondrial (JAC1)

Product Code	CSB-EP345660SVG
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P53193
Product Type	Recombinant Protein
Immunogen Species	<i>Saccharomyces cerevisiae</i> (strain ATCC 204508 / S288c) (Baker's yeast)
Purity	>85% (SDS-PAGE)
Sequence	TSTFYELFPK TFPKKLPIWT IDQSRLRKEY RQLQAQHHPD MAQQGSEQSS TLNQAYHTLK DPLRRSQYML KLLRNIDLTQ EQTSNEVTTS DPQLLLKVLVD IHDELSQMDD EAGVKLLEKQ NKERIQDIEA QLGQCYNDDK YAAAVKLTVE LKYWYNLAKA FKDWAPGKQL EMNH
Source	E.coli
Target Names	JAC1
Protein Names	Recommended name: J-type co-chaperone JAC1, mitochondrial Alternative name(s): J-type accessory chaperone 1
Expression Region	11-184
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.