



# Recombinant *Saccharomyces cerevisiae* 26S protease subunit RPT4 (RPT4)

<b>Product Code</b>	CSB-EP345691SVG-B
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	P53549
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	<i>Saccharomyces cerevisiae</i> (strain ATCC 204508 / S288c) (Baker's yeast)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	<p>           SEEQDPLLA GLGETSGDNH TQQSHEQQPE QPQETEEHHE EEPSRVDPEQ            EAHNKALNQF KRKLLRHRY DDQLKQRRQN IRDLEKLYDK TENDIKALQS            IGQLIGEVMK ELSEEKYIVK ASSGPRYIVG VRNSVDRSKL KKGVRVTLDI            TTLTIMRILP RETDPLVYNM TSFEQGEITF DGIGGLTEQI RELREVIELP            LKNPEIFQRV GIKPPKGVLL YGPPGTGKTL LAKAVAATIG ANFIFSPASG            IVDKYIGESA RIIREMFAYA KEHEPCIIFM DEVDAIGRR FSEGTSADRE            IQRTLMELLT QMDGFDNLGQ TKIIMATNRP DTLDPALLRP GRLDRKVEIP            LPNEAGRLEI FKIHAKVKK TGEFDFAAV KMSDGFNGAD IRNCATEAGF            FAIRDDRDHI NPDDLKAVR KVAEVKKLEG TIEYQKL         </p>
<b>Source</b>	E.coli
<b>Target Names</b>	RPT4
<b>Protein Names</b>	Recommended name: 26S protease subunit RPT4 Alternative name(s): 26S protease subunit SUG2 Proteasomal cap subunit
<b>Expression Region</b>	2-437
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full Length of Mature Protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.