



Recombinant *Saccharomyces cerevisiae* Uncharacterized isomerase YNR071C (YNR071C)

Product Code	CSB-MP345704SVG
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P53757
Product Type	Recombinant Protein
Immunogen Species	<i>Saccharomyces cerevisiae</i> (strain ATCC 204508 / S288c) (Baker's yeast)
Purity	>85% (SDS-PAGE)
Sequence	MSNSNGDNKY GVITIGDEKK FQATIAPLGA TLVDLKVNGQ SVVQGYSNVQ DYLTDGNMMG ATVGRYANRI AKGVFSLDDG PHKLTVNCG NTNHSSISL NLKQYKASPV ENPSKGVYVV EFKLLDDHTQ PNPNEFPGDL EVTVKYTLNV AEMTLDMEYQ AQLVRGDATP INMTNHSYFN LNKVKSEKSI RGTEVKVCSN KSLEVTEGAL LPTGKIIERN IATFDSTKPT VLHEDTPVFD CTFIIDANKD LKT TDSVSVN KLVPVFKAYH PESHKFEVS TTEPTVHLYT GDNLCGKFVP RSGFAVQQGR YVDAINRDEW RGCVLLKRGE VYTSKTQYKF DI
Source	Mammalian cell
Target Names	YNR071C
Protein Names	Recommended name: Uncharacterized isomerase YNR071C EC= 5.-.-.-
Expression Region	1-342
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.