



Recombinant *Drosophila erecta* Trypsin eta (etaTry)

Product Code	CSB-MP345791DLK
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P54629
Product Type	Recombinant Protein
Immunogen Species	<i>Drosophila erecta</i> (Fruit fly)
Purity	>85% (SDS-PAGE)
Sequence	IVG GADTTNYHTK YVVQLRRRSS PSSSYAQTCG GCILDAVTIA TAAHCVYNRE AENFLVVAGD DSRGGMSGVV VRVSKLIPHE LYNATIMDND IALVIVDPPL PLASSSTMEA IEIAAEQPAV GVQATISGWG YTKENGLSSD QLQQVNPVV DSEKQAEAYY WRPISEGMLC AGLSEGGKDA CQGDSGGPLV VANKLAGIVS WGEGCARPNY PGVYANVAYF KDWIASRV
Source	Mammalian cell
Target Names	etaTry
Protein Names	Recommended name: Trypsin eta EC= 3.4.21.4
Expression Region	28-258
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.