



Recombinant Solanum tuberosum 26S protease regulatory subunit 6B homolog

Product Code	CSB-MP345806FIG
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P54778
Product Type	Recombinant Protein
Immunogen Species	Solanum tuberosum (Potato)
Purity	>85% (SDS-PAGE)
Sequence	MATAMVLDPK PAEKL PATRP ETSITDVPSD GEDDLYARLK SLQRQLEFIE IQEEYVKDEL KNLRRHLRA QEEVKRIQSV PLVIGQFMEM IDQNNNAIVGS TTGSNYYVRI LSTINRELLK PSASVGLDRH SNALVDVLPP EADSSISLLS QSEKPDVTYN DIGGCDIQKQ EIREAVELPL THELYKQIG IDPPRGVLLY GPPGTGKTML AKAVAHHTTA AFIRVVGSEF VQKYLGEGR MVRDVFR LAK ENAPAIIFID EVDAIATARF DAQTGADREV QRILMELLNQ MDGFDQTVNV KVIMATNRAD TLDPALLRPG RLDRKIEFPL PDRRQKRLVF QVCTAKMNLG DEVDLEDYVS RPKISAAEI TAICQEAGMH AVRKNRYVIL PKDFEKG YRT NVKKPDTDFE FYK
Source	Mammalian cell
Protein Names	Recommended name: 26S protease regulatory subunit 6B homolog
Expression Region	1-413
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.