



Recombinant Escherichia coli Uncharacterized sugar kinase YgcE (ygcE)

Product Code	CSB-MP345834ENV
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P55138
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli (strain K12)
Purity	>85% (SDS-PAGE)
Sequence	MSKKYIIGID GGSQSTKVVM YDLEGNVCE GKGLLQPMHT PDADTAEHPD DDLWASLCFA GHDLMSQFAG NKEDIVGIGL GSIRCCRALL KADGTPAAPL ISWQDARVTR PYEHTNPDVA YVTSFSGYLT HRLTGEFKDN IANYFGQWPV DYKSWAWSED AAVMDKFNIP RHMLFDVQMP GTVLGHITPQ AALATHFPAG LPVVCTTSDK PVEALGAGLL DDETAVISLG TYIALMMNGK ALPKDPVAYW PIMSSIPQTL LYEGYGIRKG MWTVSWLRDM LGESLIQDAR AQDLSPEDLL NKKASCVPPG CNGLMTVLDW LTNPWEPYKR GIMIGFDSSM DYAWIYRSIL ESVALTLKNN YDNMCNEMNH FAKHVIITGG GSNSDLFMQI FADVFNLPAR RNAINGCASL GAINTAVGL GLYPDYATAV DNMVRVKDIF IPIESNAKRY DAMNKGIFKD LTKHTDVILK KSYEVMHGEL GNVDSIQSWS NA
Source	Mammalian cell
Target Names	ygcE
Protein Names	Recommended name: Uncharacterized sugar kinase YgcE EC= 2.7.1.-
Expression Region	1-492
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.